

Photron FASTCAM Analysis (PFA)

Ver.1.2 Release Notes

Ver.1.2.2.1

Release Date: 11th May, 2015

Fixed issue

- ✓ AVI files are not opened.

Ver.1.2.2

Release Date: 17th April, 2015

Fixed issue

- ✓ Under certain conditions, image files which is saved as 12bits mode are not opened.

Changes

- ✓ Bayer files which are generated from FASTCAM Multi, FASTCAM Mini WX100, FASTCAM Mini WX50 become available.

Ver.1.2.1.1

Release Date: 23rd July, 2014

Fixed issue

- ✓ Under certain conditions, the order of frames is not correct when AVI files which are generated from PhotoCam SpeederV2 are opened.

Ver.1.2.1
Release Date: 26th June, 2014

Changes

- ✓ AVI files which are generated from PhotoCam SpeederV2 become available.

Ver.1.2.0
Release Date: 16th June, 2014

New features

- ✓ Standalone launch.
- ✓ New tracking algorithm.(Binarized center of gravity)
- ✓ Tracking point can be adjusted at the zoom window.
- ✓ Scale is added on over view graph area.
- ✓ Jump to frame from a graph point.
- ✓ Selectable output files.

Changes

- ✓ Windows 8.1 is supported.
- ✓ Bayer files which are generated from FASTCAM Mini UX50 become available.
- ✓ MRAW, RAW, RAWW files which are generated from FASTCAM Mini UX50 become available.
- ✓ Tracking point is changed to a cross mark from a circle mark.
- ✓ Zoom window is fitted to template size.
- ✓ Zoom window refresh during tracking is revised to smooth.
- ✓ Time unit on the graph is changed from “ms” to “s”.
- ✓ No license PC’s functional limitation is changed.
- ✓ Language list is revised.